Material

Schaft:
Polyamid (PA)
Polyester (PES)
Polyurethan (PUR)
Zellkautschuk

Material

Upper:
Polyamide (PA)
Polyester (PES)
Polyurethane (PUR)
Cellular rubber

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Polyester (PES)
Polyurethan (PUR)

Upper- and foot plate lining:

Polyester (PES)
Polyurethane (PUR)

Sohle:

Polyurethan (PUR)

Sole:

Polyurethane (PUR)

Schutzbügel:

Stahl

Protective frame:

Steel

AUSTRIA
 Bauerfeind Ges.m.b.H.
 Hainburger Straße 33
 1030 Wien
 P +43 (0) 800 4430-130
 F +43 (0) 800 4430-131
 E info@bauerfeind.at

ITALY
 Bauerfeind AG
 Filiale Italiana
 Piazza Don Enrico Mapelli 75
 20099 Sesto San Giovanni (MI)
 P +39 02 8977 6310
 F +39 02 8977 5900
 E info@bauerfeind.it

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english

Dear Customer,

thank you for choosing a product from Bauerfeind. **GloboPed Postop** is a high-quality medical product. Please read through these **instructions for use** carefully and if you have any questions contact your doctor, medical retailer or our Technical Service department.

Product Description

The **GloboPed Postop** forefoot-relieving orthosis is a prescribable product that may only be worn under a physician's direction and only in accordance with these instructions for use and for the applications named.

Purpose

The **GloboPed Postop** forefoot-relieving orthosis is a medical product for targeted relief around the forefoot and front metatarsal area.

Indications

- For post-operative or post-traumatic conditions of the forefoot and front metatarsal area
- For diabetic/neuropathic plantar ulcers in the forefoot area
- For rheumatic feet and algodystrophic syndromes

Contraindications / Application Risks

Mobility and balance may be compromised based on the purpose and desired function of the forefoot-relieving orthosis.

To date, there are no known side effects. Correct usage and fitting by a trained specialist¹ is assumed. To date, there are no known hypersensitivities that cause illnesses. This product should be fitted and worn only after consulting your physician in the following clinical situations:

- Skin diseases / injuries in the body part being treated, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up.
- Acroparesthesia in the body parts being treated
- Compromised lymph system, circulation problems, unclear soft tissue inflammation around the product's area of application
- Raised scars with swelling, redness, and excessive heat build-up

Absolute contraindications are not known. Relative contraindications:

- Need for complete relief without touching the ground
- Need for continued mechanical stabilization of the foot or ankle
- Desintegration of the skeleton of the foot as part of a diabetic neuropathy of the hindfoot
- Limited movement of the ankle with reduced dorsiflexion

The **GloboPed Postop** forefoot-relieving orthosis may only be worn with a correctly installed protective frame that protects the forefoot from mechanical influence.

⚠ Any open wounds on the foot must be protected by a bandage or similar device. Do not ride a bicycle or drive a car or any other vehicle while wearing the orthosis. Please see your physician immediately should you notice any irregular changes in your body (i.e. increase in symptoms).

Application Instructions

Product liability cannot be assumed in the event of inappropriate use. The forefoot-relieving orthosis should be fitted while the patient is sitting and can be worn both on the left and right foot. Secure the Velcro fastenings so it fits snugly against the foot, although not too tight in order to avoid pressure points on the skin or damage to arteries and nerves. Always wear a non-skid shoe with the same heel height on your healthy foot. If necessary, a physician may prescribe that the height of the heel be changed by an orthotist. Be sure to have your **GloboPed Postop** forefoot-relieving orthosis checked regularly for wear and tear and to ensure it is working properly by an orthotist. If necessary, exchange the device. This product has been tested by our in-house quality management system. However, should you have any complaint, please contact our local expert adviser or our customer service department. Use crutches when wearing a forefoot-relieving orthosis.

Wearing the **GloboPed Postop** forefoot-relieving orthosis is only recommended in hospitals and clinics or within your own home but not in public areas. Combining the aid with other medical products (e.g. compression stockings, supports) may only be considered after having sought advice from a qualified specialist¹.

Cleaning Instructions

We recommend washing the orthosis at 30 °C by hand using a delicate detergent. We recommend closing the Velcro fastenings before washing in order to preserve their functionality and avoid damage. The orthosis should not be exposed to direct heat when drying it (i.e. heating vents, sunlight, etc.).

Accessories

Model	Size	Article Number
Protective Frame, Black Plastic	S+M	3 662490 002
Protective Frame, Black Plastic	L+XL	3 662490 003
Protective Frame, Aluminum	S-XL	3 662490 001

Notes on reuse

This product is intended for the treatment of one individual patient and may not be reused or used by others.

Warranty

The legal provisions of the country in which the product was purchased apply. Country-specific warranty regulations that govern the relationship between the retailer and the buyer also apply where relevant. If you believe that a claim may be made under the warranty, please first contact the retailer from which you purchased the product directly.

Please do not make any changes to the product yourself. This particularly applies to individual adjustments made by specialists¹. Please follow our instructions for use and care. These are based on our experience gained over many years and ensure the functionality of our medically effective products for a long time. Only aids that function optimally can help you. Failure to follow the instructions may also restrict the warranty.

First-time application and fitting of the forefoot-relieving orthosis may only be carried out by trained specialists¹, otherwise product liability cannot be assumed.

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

You can dispose of the product in accordance with national legal requirements.

Version: 2019-07

¹ A specialist is any person who is authorized according to the state regulations that are relevant to you for fitting and instruction in the use of orthopedic aids.

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